



KOO MEDICAL



Compressor nebulizer system



COMPONENTS:

Compressor, carrying case, air tubing, nebulizer, mouthpiece, nosepiece, adult mask, pediatric mask, air filter.



- Easy to use
- Easy to clean
- Easy to carry
- Complete with accessories
- Quiet during use

KOO EUROPE srl Viale delle Industrie ,5 20020 Arese (MI) Italy
Tel. 0039 02 93589495 info@kooeurope.com - www.koomedical.com



KOO MEDICAL



Compressor nebulizer system



Nebulizer



Mouthpiece



Nosepiece



Air tubing



Adult mask,
pediatric mask



Carrying case



Air filter

Voltage	230V/50Hz – 120V/60Hz
Power consumption	180 VA
Nebulizer capacity	5 ml
Particles size	0,5 to 10 μ m
MMAD	4 μ m
Nebulization rate	\geq 0,2ml/min.
Noise level	60 dBA a 1 mt
Compressor's pressure	\geq 29 PSI (\geq 0,20 Mpa)
Operating pressure	8-17Psi (0,056 – 0,117 Mpa)
Air flow	4~8 LPM
Operating temp.	10°C / 40°C (50°F/104°F)
Storage Temperature:	-20°C max. 70°C (-4°F/158°F)
Humidity:	10% - 95%
Dimensions:	203x153x115 mm
Weight:	1,4 kg

KOO EUROPE srl Viale delle Industrie ,5 20020 Arese (MI) Italy
Tel. 0039 02 93589495 info@kooeurope.com - www.koomedical.com