

CRYSTAL AIRWAY MASK™



- **Latex - Free**
- **Available in 7 sizes**
- **Different Colors for Sizes**
- **PVC Material**
- **Disposable**

Features and Benefits

- **Rotationally-molded laryngeal cuff**

The process of rotational molding produces a cuff with a stronger yet thinner wall as compared to other methods of cuff manufacture. This process yields a cuff that is ultra-soft next to the delicate tissues of the larynx and the laryngeal aditus, reducing the potential for trauma and creating a more custom seal. This custom seal capability allows greater ventilation options and more efficient airway protection without the dangers of high cuff pressures.

- **Reinforced posterior aspect of cuff**

To aid in insertion of the Crystal Airway Mask™ the posterior aspect of the cuff (away from the delicate laryngeal tissues) is reinforced with a thicker wall than the patient side of the cuff. This reinforcement reduces the possibility of the distal portion of the cuff folding back on itself during insertion. This type of reinforcement is possible only with rotational molding, as it is produced in entirely one contemporaneous piece.

- **Integrated inflation lumen**

Unlike other devices on the market, the Crystal Airway Mask™ cuff inflation line is integrated into the wall of the ventilation tube. This serves to protect the inflation line and reduce the possibility of leaks caused by erosion against or puncture by the teeth.

- **Single use device**

The Crystal Airway Mask™ is a single use device, intended for disposal after use. Today's environment of resistant pathogens and their role in nosocomial infections demand disposability of all products that come into contact with mucous membranes.

CRYSTAL AIRWAY MASK™

Technical Characteristics

Product Number	Mask Size	Patient Weight (Kg) Category	I.D. Tube (mm)	O.D. Tube (mm)	Maximum Cuff Volume (ml)	Quantity Per Box
KM-877	1	<5	5.6	8.6	5	10
KM-867	1.5	5 - 10	6.4	10	7	10
KM-857	2	10 - 20	7.8	12	10	10
KM-847	2.5	20 - 30	9.5	14	15	10
KM-837	3	30 - 50 Small Adult/Children	10	16	20	10
KM-827	4	50 - 70 Adult Female/Male	11.4	18	35	10
KM-817	5	>70 Adult Male	11.9	20	48	10



Color - coded

Distributed by:



KOO MEDICAL EQUIPMENT

KOO EUROPE SRL

Via della Praia 12/D, 10090 - Buttigliera Alta (TORINO), ITALY
Tel : +39 011 934 8192 Fax : +39 011 934 8627
info@kooeurope.com

KOO MEDICAL JAPAN

3-22-11-2F, Hongo, Bunkyo-ku, Tokyo, JAPAN
Tel : +81 3 3818 6635 Fax : +81 3 3818 6632

KOO ASIA LTD.

Unit 23, 9 Floor Goodwill Ind. Building, No. 36-44,
Pak Tin Pa Street, Tsuen Wan, HONG KONG
Tel : +852 2110 4780 Fax : +852 2110 4783
sales@koo-asia.com

KOO AMERICAS

sales@koo-americas.com

MKT_KM-1008