



Innovative range of masks
N.I.V. Bluestar Plus



The air at service of breathing

Lightness

Seal

Comfort

KOO MEDICAL EQUIPMENT
www.koomedical.com

The air at service of breathing...

The lightness of the air
The seal of the air cushion
The comfort of the anatomic cushion

NEW

Facial Mask - Intensive

Bi level CPAP **Non-vented**

Its design makes it widely versatile, the presence in the Kit of 3 connectors allows it to operate with different configurations of circuits. The connector with hole for exhalation allows to transform the mask from **Non-vented** to **vented**. A port on the connector allows the O₂ enrichment or monitoring.



BLUESTAR PLUS Intensive

- *Replaceable anatomic air cushion*
- *Wide range of fixing points*
- *Headgear mask fixer in highly breathable material and with wide bearing surface for better forces distribution*
- *Operative openings for O₂ enrichment and for monitoring*
- *Port for the passage of the gastric tube, or the introduction of the bronchoscope, or for simple needs of the patient*
- *Front support fully adjustable and removable*



The air at service of breathing...

*The lightness of the air
The seal of the air cushion
The comfort of the anatomic cushion*

NEW

Facial Mask - Secure

Bi level CPAP **Non-vented**

The presence of a safety valve Anti-**Asphyxia**, permits a reasonable airflow even in case of the flow lack coming from the fan.

There is a rotating connector a 360° blue color with a diameter of 22 mm F.



BLUESTAR PLUS Secure with Security valve A.A.V.

- *Replaceable anatomic air cushion*
- *Wide range of fixing points*
- *Headgear mask fixer in highly breathable material and with wide bearing surface for better forces distribution*
- *Operative openings for O₂ enrichment and for monitoring.*
- *Port for the passage of the gastric tube, or the introduction of the bronchoscope, or for simple needs of the patient*
- *Safety valve with immediate opening in the event of flow lack*
- *Front support fully adjustable and removable*

The air at service of breathing...

The lightness of the air
The seal of the air cushion
The comfort of the anatomic cushion

NEW

Facial Mask - Home CARE

Bi-level CPAP **Vented**

The presence of exhalation holes, the safety valve Anti Asphyxia and the replaceable air cushion, make this mask the most versatile for use in the ventilation of the patient at home-care.



BLUESTAR PLUS Home-CARE with exhalation holes and security valve A.A.V Home-Care

- *Replaceable anatomic air cushion*
- *Wide range of fixing points*
- *Headgear mask fixer in highly breathable material and with wide bearing surface for better forces distribution*
- *Operative openings for O₂ enrichment and for monitoring.*
- *Port for the passage of the gastric tube, or the introduction of the bronchoscope, or for simple needs of the patient*
- *Safety valve with immediate opening in the event of flow lack*
- *Exhalation holes for a correct evacuation of the exhaled breath*
- *Front support fully adjustable and removable*

The air at service of breathing...

Our range



Bluestar Plus Intensive

KM - 310 B	Extra Large	Each kit includes: Mask with air cushion, transparent connector 22 mm M, blue connector 22 mm F, exhalation valve, headgear mask fixer, bulb pump.
KM - 311 B	Large	
KM - 312 B	Medium	
KM - 313 B	Small	
KM - 314 B	Pediatric	



Bluestar Plus Secure AVV

KM - 310 V	Extra Large	Each kit includes: Mask with air cushion, blue connector 22 mm F with security valve, headgear mask fixer, bulb pump.
KM - 311 V	Large	
KM - 312 V	Medium	
KM - 313 V	Small	
KM - 314 V	Pediatric	



Bluestar Plus Home CARE

KM - 310 E	Extra Large	Each kit includes: Mask with air cushion, transparent connector 22 mm M with security valve and exhalation holes, headgear mask fixer, bulb pump.
KM - 311 E	Large	
KM - 312 E	Medium	
KM - 313 E	Small	
KM - 314 E	Pediatric	



Accessories

KM - 3101	Extra Large	Cushion replacement
KM - 3111	Large	
KM - 3121	Medium	
KM - 3131	Small	
KM - 3141	Pediatric	
KM - 3104	Transparent connector 22 mm M	
KM - 3114	Blue connector 22 mm F	
KM - 3113	Forehead Space	
KM - P009	Headgear Large	
KM - P010	Headgear Small	
KM - 903	Exhalation valve	



The air at service of breathing...

The lightness of the air

The use of the air cushion allows a reduction of weight.
A mask lighter requires a lower pulling force
from the headgear, to keep it stable

The seal of the air cushion

The air cushion, anatomically designed, ensures
an excellent control of the seal,
avoiding losses peripherals unpleasant

The comfort of the anatomical cushion

The anatomical air cushion allows a distribution
of the pressure sealing eliminating dangerous pressure's drops.
The complete adjustment of the mask and
the use of a headgear made of breathable material,
makes the mask highly accepted by most patients,
preventing the interruption of ventilation, during use of the mask.

**The best handy and low-cost
solution for non-invasive
ventilation with a device that
allows the replacement of
the cushion and the shell**



Responsible for **EMEA (Europe Middle East Africa) Market**
Koo Europe s.r.l.

viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market**
Koo Asia Ltd

Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong
Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Koo Medical Equipment (Shanghai) Co.,Ltd.
www.koomedical.com

1 Songjiang Industrial Park - Songjiang, Shanghai, 201613 P.R.China