



www.koomedical.com

Aerosol **FISIOJET** kit with:



- High quality compressor
- Nebulizer Fisiojet KM 1000 with high performance in polycarbonate
- Latex free DHEP and BPA free
- Elbow connector
- Mouthpiece
- Nose piece
- Connection tube
- In a carry bag

Every item is in a box

Code:

FJ001 Junior with Fisiojet infant Mask and Fisiojet pediatric Mask

FJ002 Senior with Fisiojet Adult Mask



FISIOJET nebulizer for aerosol therapy

- High efficiency nebulizer with inspiratory and expiratory valves
- This type of supply optimizes the inhalation therapy, guaranteeing, maximum medication absorption during inspiration and minimum medication loss during expiration
- Minimal flow 3,0 l/min maximum flow 6,0 l/m
- Minimum cap. 2 ml maximum cap. 8 ml
- In polycarbonate DHEP and BPA free
- Latex free

Provided with the following items:
elbow connector, mouthpiece, nose piece

- single patient use

Code:

KM-1002 in a box



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel: +852 2110 4780 Fax: +852 2110 4783 - Email: sales@koo-asia.com


FISIOJET JUNIOR

nebulizer for aerosol therapy with:

- Nebulizer, elbow connector, mouthpiece, nose piece
- Fisiojet infant mask, Fisiojet pediatric mask
- Connecting tube ca. 120 cm
- Single patient use

Code:

KM-1003 in a box





FISIOJET SENIOR

nebulizer for aerosol therapy with:

- Nebulizer, elbow connector, mouthpiece, nose piece
- Fisiojet adult mask
- Connecting tube ca. 120 cm
- Single patient use

Code:

KM-1004 in a box



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com
Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

FISIOJET

facial mask for aerosol therapy



- In PVC DEHP and BPA free
- Latex free
- Soft and comfortable
- Minimum drug loss
- A valve for the exhalation of CO₂
- Single packaging
- Single patient use

code	size	items per box
KM-1005	adult	50
KM-1006	small pediatric	50
KM-1007	large pediatric	50
KM-1008	infant	50

Accessories:

code	items per box
KM-1009 Fisiomet nebulizer	50
KM-1010 standard elbow connector	50
KM-1011 nose piece to use with the mouthpiece	50
KM-1012 standard mouthpiece	50
Km-1013 tube of aerosol ca. 120 cm	100

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23.9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com